

NMR_Service: Tipps&Tricks


I. Basics about sensitivity in NMR Spectroscopy

- The signal-to noise (S/N) of a FT NMR spectrum does only increase with the **square-root of the number of scans (and thus the time) !**


$$\text{NMR Spectroscopy: } S/N \sim \sqrt{t}$$

- So, if we have **2 times less substance**, an identical spectrum (typically a ^{13}C -spectrum) takes **4 times as long!**
- Corollary: A **2-fold increase in S/N** requires **4 times as much time!**

NS = 256



NS = 1024



400 MHz ^{13}C spectrum of a 50 mmol sample of Cyclosporine. Bottom: 1024 scans, experiment time: 40 min. Top: 256 scans, experiment time: 10 min. Clearly, the signal-to noise (S/N) only increases with the **square-root of the number of scans (and thus the time) !**

II. How much substance do I need?


A. Open-Access 300 MHz: characterization incl. ^{13}C -spectrum and/or DEPT135, APT

Since the length of ^{13}C -based experiments is limited to **40 minutes** on both 300 MHz spectrometers, you need **0.6 mL** of a **50 mM solution** (at least) of your molecule to get a good spectrum.

For a typical organic molecule this is about **3-15 mg**:

Molecular mass (g)	100	200	300	400	500
Required quantity (mg) for 0.6 mL for a 100 mM solution	3	6	9	12	15

An example of what we consider a good ^{13}C spectrum is shown below for a 50 mmol cyclosporine sample.



300 MHz ^{13}C spectrum of a 50 mmol sample of Cyclosporine. 1024 scans, experiment time: 40 min.

B. Service: Complete characterization incl. ^{13}C -spectrum and a set of 2D spectra

1. ^{13}C -spectrum


A very important fact: although it may seem counter-intuitive, the ^{13}C spectrum largely determines the length of the NMR-characterization! This characterization will be even longer if a DEPT and/or APT is also requested!

The reason is that the 2D spectra can be kept short because the large signal-to-noise on the ^1H channel allows us to use the minimum number of scans, even for relatively low sample concentrations!

Note that very often the ^{13}C -spectrum is not even needed for complete characterization because one can read off the ^{13}C -shifts from the indirect dimension of the HSQC and HMBC spectra.

Important: If you still want a 1D ^{13}C spectrum, we strongly recommend the APT (DEPTQ) spectrum (see below), as it contains everything you need for analyzing ^{13}C spectra.

^{13}C + DEPT or ^{13}C + APT or the completely needless combination ^{13}C + DEPT + APT are just redundant and unnecessarily lengthen the overall measurement time (and increase the bill for your boss...).



400 MHz ^{13}C -APT spectrum of Strychnine. The signals of CH and CH_3 are *negative*, while CH_2 and quaternary carbons including the solvent carbon are *positive*.

2. Examples

As an example, the following settings provide good spectra (see below) for a **50 mM sample** on one of our 400 MHz spectrometers:


- ^1H : 8 scans (1')
- ^{13}C -APT (DEPTQ): 1024 scans (40')
- COSY, 128 increments: 2 scans per increment (4,5')
- HSQC, 128 increments: 2 scans per increment (4,5')
- HMBC, 128 increments: 4 scans per increment (10')

Total experiment time is **60 minutes**

66% of the measurement time is spent for the ^{13}C or APT spectrum...


➤ This percentage increases with decreasing sample concentration...

50 mmol Cyclosporine, C6D6
NS = 1024
40 min




400 MHz ^{13}C spectrum of a 50 mmol sample of Cyclosporine. 1024 scans, experiment time: 40 min.

NS = 2




400 MHz 2D ^1H - ^{13}C edited HSQC spectrum of a 50 mmol sample of Cyclosporine. 128 t_1 points, 2 scans, experiment time: 4,5 min. The signals of CH and CH₃ are positive (black

contours), while CH_2 are negative (red contours). Quaternary carbons do not provide signals in HSQC spectra!



Zoom of the 400 MHz 2D 1H - ^{13}C edited HSQC spectrum of a 50 mmol sample of Cyclosporine. All ^{13}C -shifts can be very easily read off the from the indirect dimension.

NS = 4



400 MHz 2D ^1H - ^{13}C HMBC spectrum of a 50 mmol sample of Cyclosporine. 128 t_1 points, 4 scans, experiment time: 10 min. *Quaternary carbons* do provide signals in HMBC spectra!

C. I have only enough for a 10 mM sample. Can I still get decent data?

Yes, but the necessary measurement time will be significantly longer (and so the costs for your boss...), especially if you absolutely want the ^{13}C APT(DEPTQ) spectrum.

Remember that compared to the 50 mmol sample, the APT spectrum now needs to be recorded with 25'600 scans to achieve a comparable result !


➤ **This takes approximately 17 hours !.**

If you (and your boss...) can make do with a lower-quality ^{13}C -APT spectrum (more noise), a ^{13}C -APT spectrum on a 10 mM sample can of course be recorded in 6 hours. The resulting signal-to-noise (intensity of the resonance, basically) will be ~half as much as compared to a 17-hour spectrum.

An example of what we consider a decent ^{13}C spectrum is shown below for a 9 mmol cyclosporine sample.

- decent: all carbons are visible, but the S/N is sometimes poor, and the noise is intense.

8.9 mmol Cyclosporine
NS = 4096
3h 20 min



400 MHz ^{13}C spectrum of a 9 mmol sample of Cyclosporine. 4096 scans, experiment time: 3h20 min.

At 400 MHz, if you expect a decent ^{13}C -spectrum, we recommend a concentration of no less than 10 mM of your molecule in a deuterated solvent.

For a typical organic molecule this is about **0.6-3 mg**.

Molecular mass (g)	100	200	300	400	500
Required quantity (mg) for 0.6 mL for a 10 mM solution	0.6	1.2	1.8	2.4	3

All proton-detected experiments including COSY, HSQC and HMBC are still possible at relatively low concentrations!

For a 10 mM sample we recommend ns=4, ns=8 and ns=16 for COSY, HSQC and HMBC, respectively.

For a 5 mM sample we recommend ns=8, ns=16 and ns=32 for COSY, HSQC and HMBC, respectively.

At 2.5 mM a good HMBC will take too long, but COSY and HSQC will still be of good quality.

**D. I have only a 5 mmol solution of my compound.
Can I get a ^{13}C spectrum? 2D spectra?**

- For the ^{13}C spectrum, if you consider what we have said above, the answer is no, unless:
You're willing to get a (very)poor carbon spectrum.
- For the 2D spectra, the answer is yes:
It means that you should be prepared to characterize your compound only with 2D spectra (COSY, HSQC, HMBC with limitations). Note again that most of the time the ^{13}C -spectrum is superfluous and not even needed for complete characterization because one can read off the ^{13}C -shifts from the indirect dimension of the HSQC and HMBC spectra.